

# Speech in Noise Testing

## An expanded and integrated view

Francis Kuk, Ph.D.  
ORCA-USA, WS Audiology

# SPEAKER DISCLOSURE

**Relevant Financial Relationships: Employee of WS Audiology**

**Relevant Non-financial Relationships: None**



# AGENDA



- Realistic communication model
  - Objective vs subjective intelligibility assessment
  - Working memory limits speech in noise comprehension
  - Role of cognitive factors , e.g. listening effort
  - Psychological factors to stay in noise
    - noise acceptance, including the TNT
- Repeat-Recall Test (RRT) and Quick-RRT
- Clinical use for the Quick RRT



# REALISTIC COMMUNICATION

## Listener activities

We listen to  
understand



## Clinicians' task

Need to understand how much they  
**understand** and how much they  
**think they understand**

# OBJECTIVE VS SUBJECTIVE INTELLIGIBILITY TESTS

## Objective tests (actual understanding)

- Listeners repeat sentences/ words they hear or understand
- Takes some time
- All clinical speech tests we use (NU-6, HINT, QuickSin) – single nonsense words to meaningful sentences
- Somewhat predictable from audibility (may have cognitive component)

## Subjective tests ("think" understanding)

- Listeners estimate how much they think they heard or understand
- Takes less time (in principle)
- HINT and QuickSin have been adapted to yield subjective SRT50 – sentences
- Typically, we expect subjective to be similar to objective (for words and sentences).

# SCORING THE QUICKSIN AS A SUBJECTIVE INTELLIGIBILITY TEST (OU & WETMORE 2020)

List 1		Score
1. A <u>white</u> <u>silk</u> jacket goes with <u>any</u> shoes.	S/N 25	_____
2. The <u>child</u> <u>crawled</u> into the <u>dense</u> grass.	S/N 20	_____
3. <u>Footprints</u> <u>showed</u> the <u>path</u> he <u>took</u> up the <u>beach</u> .	S/N 15	_____
4. A <u>vent</u> near the <u>edge</u> brought in <u>fresh</u> air.	S/N 10	_____
5. It is a <u>band</u> of <u>steel</u> <u>three</u> <u>inches</u> <u>wide</u> .	S/N 5	_____
6. The <u>weight</u> of the <u>package</u> was <u>seen</u> on the <u>high</u> scale.	S/N 0	_____
27.5 - TOTAL = _____ SNR Loss	<b>TOTAL</b>	_____

## OBJECTIVE TEST

Repeat the sentence and score the number of correct words. Add up total score.

## SUBJECTIVE TEST


Estimate how much of the sentence is understood using a number from 0 to 5 with 0 means nothing and 5 every word of the sentence. Add up total score.



# WHY IS IT IMPORTANT TO KNOW WHAT LISTENER THINKS S/HE UNDERSTANDS?

Perception can be affected by emotions, cognition in addition to acoustics

Perception can affect decisions (such as HA satisfaction and purchase)



Subjective  
intelligibility




Objective  
intelligibility



# PROBLEMS WITH SUBJECTIVE INTELLIGIBILITY NOT BEING THE SAME AS OBJECTIVE INTELLIGIBILITY

- Subjective = objective – reasonable expectations, expected outcome from intervention ( $S = O$ )
- Over-estimator – thinks s/he hears better than s/he actually does – denies problem, less willing to accept help ( $S > O$ )
- Under-estimator – thinks s/he has more problems than the hearing loss or speech test results suggest – critical of hearing (or intervention), difficult to please, lower satisfaction ( $O > S$ )


# SUBJECTIVE-OBJECTIVE INTELLIGIBILITY DIFFERENCE (SRT50) BETWEEN NORMAL-HEARING AND HEARING-IMPAIRED (UNAIDED MODE)



For Normal-Hearing listeners, average subjective  $SRT_{50}$  is same as objective  $SRT_{50}$  around -6 dB

For Hearing-Impaired listeners, average subjective  $SRT_{50}$  is better (or lower) than objective  $SRT_{50}$  by 1 – 1.5 dB (average around -2.5 to -3.0 dB)

# SUBJECTIVE-OBJECTIVE INTELLIGIBILITY DIFFERENCE (SRT50) BETWEEN UNAIDED AND AIDED MODES IN HEARING-IMPAIRED






For the Hearing-Impaired listeners,

- Unaided – subjective SRT > objective SRT about 1.5 dB
- Aided – subjective SRT = objective SRT (+/-0.5-1, not significant); SRT about -3 to -4 dB

# OUR RECENT STUDY ON SUBJECTIVE VS OBJECTIVE INTELLIGIBILITY SHOWS

- We used related sentences created from our TNT test passages
- Normal hearing listeners showed similar performance between subjective and objective tests
- Hearing impaired listeners showed better subjective performance than objective performance in unaided mode (almost half, based on criterion from NH)
- Hearing impaired listeners showed less difference between subjective performance and objective performance in aided mode (similar to normal hearing)

# INCIDENCE OF OVER- & UNDER-ESTIMATOR IN HI LISTENER RE:NH LISTENER



# IMPLICATIONS AND CHALLENGES

- Implications
  - Unaided - HI perceives less difficulty than reality ( $S > O$ )
  - Aided – HI's subjective intelligibility closer to NH
    - HAs restore some level of “normal” perception of intelligibility or reality
    - HI perceives less benefit than reality (benefit = aided – unaided performance)
- Challenges
  - Convincing HI listeners to try amplification when they perceive less difficulty than they really have and to show them more benefit than they perceive - **demonstration**
  - Include subjective evaluation – **as a supplement to objective measure**

# REALISTIC COMMUNICATION

## Listener activities

We listen to  
understand




## Clinicians' task


Need to understand how much they  
**understand** and how much they  
**think they understand**



Need to understand size of **working  
memory**



# WHAT IS WORKING MEMORY?

- Memory that **TEMPORARILY STORES** information for **SUBSEQUENT PROCESSING**
  - Examples
    - Following instructions
    - Map directions
  - Can be **verbal or nonverbal** (may or may not be same size)
  - Age effect, not same as cognitive impairment
  - Please visit [www.orca-us.info](http://www.orca-us.info) for “**RRT 1: Rationale and Development**” to hear a distinction between dementia, MCI and working memory
- 





# PEOPLE WITH BETTER vs POORER WORKING MEMORY



- Better working memory
  - **Better** speech in noise, lower listening effort, more willing to spend more time in noise
- Poorer working memory
  - Greater need for noise management strategies
  - Greater need to preserve signal naturalness (use of slow compression, low delay etc)



# REALISTIC COMMUNICATION

## Listener activities

We listen to understand



We try to retain what we can understand



We integrate information & formulate response



## Clinicians' task

Need to understand how much they **understand** and how much they **think they understand**

Need to understand size of **working memory**

Need to measure **perception** of **effort** and **willingness to stay in noise**




# WHAT IS LISTENING EFFORT?

- “...*the deliberate allocation of mental resources to overcome the obstacle in goal pursuit when carrying out a listening task...*” (Pichora-Fuller et al, 2016)
- Means to measure
  - Physiological
    - Heart Rate Variability
    - Pupillometry
    - EEG (alpha)
  - Behavioral
    - Dual task
    - Rating




LISTENING  
EFFORT CAN BE  
*BIASED* BY  
PERFORMANCE;  
BUT NOT  
ALWAYS



A subjective task;  
As a cross check on performance;  
May be especially sensitive when performances are similar

# WHY NOISE ACCEPTANCE MEASURE (WHILE MAINTAINING SPEECH UNDERSTANDING) MAY BE BETTER AT ESTIMATING COMMUNICATION?




# A NOISE ACCEPTANCE TEST FROM ORCA

# TRACKING OF NOISE TOLERANCE (TNT)

- Measures to ensure reliability
  - Fixed speech input level
  - Tracking noise level for 2 min
  - Use multiple equivalent passages
- Measures to be more real-life
  - Speech filtered according to input level to approximate speech spectra of increased vocal effort
- Measures to ensure intelligibility as main criterion
  - Specify intelligibility criteria – > 90%
- Direct interpretation
  - **TNT** = TNL – speech
  - Higher TNT, greater noise acceptance

# TRACKING OF NOISE TOLERANCE (TNT) - PHYSICAL SETUP



Headphone  
(unaided) also


# INSTRUCTIONS – TOLERABLE NOISE LEVEL

- You will hear some noise in the background while you listen to the male talker. The noise will automatically get louder. I want you to monitor the noise level and maintain the loudest noise level you can put up with while still understanding 90% of the words in the story. If the noise becomes too loud, where you can no longer put up with it or understand less than 90% of the words in the story, you can turn the noise down by pressing and holding the space bar. If it appears softer than before, you should allow the volume to increase by letting go of the space bar. If it is louder than before, you should turn the volume down to keep at the same level by pressing the space bar again. Your ability to understand speech should never change to below 90%. The test will run for two minutes and then stop.*



# TNL TRACKING RESULTS (INDIVIDUAL)

TNT test



Speech  
Noise

Instantaneous noise (dB SPL): 79  
Time remaining (sec): 0

Quick TNL (dB SPL): 77.7  
Quick TNT (dB): 2.7  
Time to stabilize (sec): 16.9  
Average TNL (dB SPL): 78.7  
Average TNT (dB): 3.7  
Average peak (dB SPL): 79.6  
Average valley (dB SPL): 77.2  
Average excursion (dB): 2.4

Print

Save

Quit


New session

Tracking of Noise  
Tolerance (TNT)

Version 0.03  
Copyright © 2023

# SUBJECTIVE SPEECH INTELLIGIBILITY DURING TNT TRACKING

TNT test



Peaks - about 30-50%  
Average - about 60-70%  
Valleys - about 90%

Speech	
Noise	
Instantaneous noise (dB SPL):	79
Time remaining (sec):	0
Quick TNL (dB SPL):	77.7
Quick TNT (dB):	2.7
Time to stabilize (sec):	16.9
Average TNL (dB SPL):	78.7
Average TNT (dB):	3.7
Average peak (dB SPL):	79.6
Average valley (dB SPL):	77.2
Average excursion (dB):	2.4

- Print
- Save
- Quit
- New session




The Tracking of Noise Tolerance (TNT) test is a test of *subjective speech intelligibility and more*

Specifically, it tells us the noise level that a listener can accept to understand different amount of speech based on his/her criterion


# PEOPLE WHO ARE **SATISFIED** WITH THEIR HEARING AIDS ACCEPT SIMILAR AMOUNT OF NOISE AS PEOPLE WITH **NORMAL HEARING**

- Seper et al (2018) found in 17 satisfied wearers of hearing aids that their TNL were:
  - Babble noise = 79.2 dB (vs normal 79.5)
  - Continuous noise = 77 dB (vs normal 76.5)
  - ALL but one subject had TNL much below 70 dB SPL




**WEARERS WHO  
ARE MORE  
SATISFIED WITH  
THEIR HEARING  
AIDS ACCEPT  
MORE NOISE**

MarkeTrak III - Overall Hearing Aid Satisfaction Rating



MarkeTrak III - Proportion of Noisy Situations Satisfied



# THE USEFULNESS OF THE TRACKING OF NOISE TOLERANCE (TNT) TEST

It offers a *new dimension* of outcome measurement (i.e., *individual criterion of subjective intelligibility*)

SNR that listener thinks s/he needs to understand > 90% of passage materials of simple content based on his/her *own criterion* of 90%


Profiling

Fine-tuning and comparison of HAs/features


Measurement/prediction of HA satisfaction/success



# REALISTIC COMMUNICATION



# WE WANT TO KNOW THESE ABOUT OUR PATIENTS




So you can profile your patient, gauge/predict his success with intervention such as hearing aids, rehab training etc, intervene appropriately and optimally, and validate performance/success of intervention





# WHAT DO YOU MEASURE IN YOUR CLINIC?

- Objective speech test in **quiet?**
    - NU-6, W-22, CCT, NST etc
  - Objective speech test in **noise?**
    - QuickSin, HINT, CST etc
  - **Subjective** speech intelligibility measurement (in noise)?
    - QuickSin, CST
  - Assessment of **working memory** (and or cognition)?
    - Reading (listening) span, MoCA, digit span etc?
  - Assessment of listening **effort?**
    - Subjective rating, pupillometry etc?
  - Assessment of **noise acceptance?**
    - ANL, TNT
- 

**TIME**

A large, thick red circle with a diagonal slash from the top-left to the bottom-right, crossing out the word 'TIME'. The word 'TIME' is written in a bold, black, sans-serif font, centered within the circle. The background is white. There are faint, repeating watermarks of the word 'featurepics' and a small icon across the image.



**THE REPEAT-RECALL  
TEST (RRT) ALLOWS YOU  
TO GATHER ALL THE  
ABOVE INFORMATION IN  
A REASONABLE AMOUNT  
OF TIME**

# LATEST REPEAT RECALL TEST (VERSION 0.4.1)

REPEAT AND RECALL TEST  
RRT

ORCA *USA*

Beta Version 0.4.1.0  
ORCA-USA, WS Audiology

RRT RRT :: Setup

Setup Import Info

Administrator info

Tester ID:  Clinic name:

Subject info

ID:  Age:  MoCA score:

Audiogram (dB HL)

	.25	.5	1	2	4	8 (kHz)
Right	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>
Left	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>

HA info:

Test setup

Noise:  Babble 0°  Babble 180°  Continuous 0°  Continuous 180°

Testing:  Standard  Quick  Custom  Practice

Quiet  SNR = 15 dB  SNR = 10 dB  SNR = 5 dB  SNR = 0 dB

Passage 5  Passage 2  Passage 7  Passage 1  Passage 3

Books & Movies

# REPEAT-RECALL TEST (RRT) - CONSTRUCTION

5 themes  
(4<sup>th</sup> grade)

FOOD AND COOKING

MOVIES AND BOOKS

SHOPPING

SPORTS

MUSIC



7 lists/passages

High context

Keep the **ice cream** in the **freezer**.  
The **chef cooks** food in a **restaurant**.  
The **barbecue grill** used **hickory wood**.  
**Wash** the **fruit** in the **sink**.  
The **tart pie** had too much **lemon**.  
He tried **new foods** in **different countries**


6 sentences

Low context

Keep the **ice foods** in the **lemon**.  
The **cream cooks** food in a **country**.  
The **barbecue chef** used **hickory freezer**.  
**Wash** the **grill** in the **restaurant**.  
The **tart fruit** had too much **wood**.  
He tried new **pie** in **different sinks**

- A Software program on Windows
- Speech front at 75 dB SPL.
- Noise options:
  - 2-talker babble & **speech shaped noise**
  - 0° and **180°**
- Each list tested at one SNR = 0, 5, 10, 15, Q
- Testing can be conducted unaided and aided in sound field

# PHYSICAL SET-UP – STANDALONE OR VIA AUDIOMETER/BOOTH



# STEP 1 – REPEAT EACH SENTENCE

## REPEAT EACH SENTENCE

Sing along with the chapter nothing.

sing

chapter

nothing

His biography appeared on the back feelings.

biography

back

feelings

The surprise shared her honest quote.

surprise

shared

honest

quote

Each portrait began with an historic audience.

portrait

historic

audience

The poet revealed cover but the musical.

poet

cover

musical

The classic screamed at the big truth.

classic

screamed

big

truth

Repeat test

In progress

Recall test

Subjective rating

Save repeat  
test results

# STEP 2 – RECALL ALL (SIX) 6 SENTENCES

## RECALL AS MANY WORDS AS YOU CAN

Sing along with the chapter nothing.

sing

chapter

nothing

His biography appeared on the back feelings.

biography

back

feelings

The surprise shared her honest quote.

shared

honest

Each portrait began with an historic audience.

historic

audience

The poet revealed cover but the musical.

cover

musical

The classic screamed at the big truth.

big

truth

Repeat test

Completed

Recall test

In progress

Subjective rating

Repeat test  
results saved

Save recall  
test results



# STEP 3: RATE LISTENING EFFORT AND REPORT TOLERABLE TIME

### EFFORT RATING AND TOLERABLE TIME

	Listening effort	
Give up	<input type="button" value="11"/>	Tolerable time (0 ... 120 minutes) <input type="text" value="2"/>
Very effortful	<input type="button" value="10"/>	
	<input checked="" type="button" value="9"/>	
	<input type="button" value="8"/>	
	<input type="button" value="7"/>	
	<input type="button" value="6"/>	
Moderately effortful	<input type="button" value="5"/>	
	<input type="button" value="4"/>	
	<input type="button" value="3"/>	
	<input type="button" value="2"/>	
Not effortful	<input type="button" value="1"/>	

Repeat test  
Completed


Recall test  
Completed

Subjective rating  
In progress

YOU HAVE  
COMPLETED  
TESTING AT  
ONE SNR ...


REPEAT THE  
PROCESS AT  
ANOTHER  
SNR

Maximum administration time  
– 5 SNRs by 2 contexts or 20-  
25 minutes

		Repeat (%)	Recall (%)	Effort (1 - 11)	Tolerable time (min)
Completed		70	30	9	2
Start		-	-	-	-
Start		-	-	-	-
Start		-	-	-	-
Start		-	-	-	-
Start		-	-	-	-
Start		-	-	-	-
Start		-	-	-	-
Start		-	-	-	-
Start		-	-	-	-


P-I Profile Print Save data New session Exit

# THE P-I FUNCTION ALLOWS ESTIMATION OF REPEAT SNR AT DIFFERENT CRITERIA




Estimating SRT at higher criteria (75%, 85%) yields more realistic appraisal of SNR needs in daily COMMUNICATION situations

# THE PATIENT'S PERFORMANCE IS COMPARED AGAINST REFERENCE DATA





Layman's terms for SNR  
(Wu et al, 2018)


# DISPLAY OF AVAILABLE INFORMATION



Context	SNR (dB)	Repeat (%)	Recall (%)	Effort (1-11)	Tolerable time (mins)
High	0	5	5	9	3
	5	25	25	9	5
	10	70	40	6	8
	15	100	60	2	20
	Quiet	100	65	2	58
Low	0	0	0	10	1
	5	25	5	9	1
	10	60	35	9	2
	15	85	20	6	5
	Quiet	95	20	6	12


Context	SRT-50 (dB)	SRT-75 (dB)	SRT-85 (dB)
High	7.8	10.4	12
Low	8.7	12.3	14.8







FeaturePics.com - I1269436


# THE QUICK REPEAT-RECALL TEST




# QUICK RRT TEST FLOW



"Repeat each sentence as you hear it."




"On a scale of 1-10, how effortful was it to hear the speech."



"How long would you spend in this environment?"




# QUICK RRT DISPLAY



- Excellent
- Very good
- Good
- OK, improvement desirable
- Poor

Context	SNR (dB)	Repeat (%)	Recall (%)	Effort (1-11)	Tolerable time (mins)
High	0				
	5				
	10				
	15	95	55	3	79
	Quiet				
Low	0				
	5	70		8	20
	10	75		7	30
	15	100		5	50
	Quiet				

Context	SRT-50 (dB)	SRT-75 (dB)	SRT-85 (dB)
High	-	-	-
Low	5	7.8	11.2



# WHAT ARE THE **DIFFERENCES** BETWEEN THE FULL-RRT AND Q-RRT?

## Full-RRT


- Takes 20-25 minutes to complete
- Dual task (SiN + recall)
- Allows study of speech in noise, recall, LE, TT and use of ***context*** at all 5 SNRs
- LE and TT rating biased by recall
- Research tool

## Quick-RRT

- Takes 6 minutes to complete
- Single task (i.e., SiN, recall)
- Allows study of speech in noise, LE, TT at 3 SNRs; recall at 1 SNR only
- LE and TT based on SiN perception
- Clinical tool


# RECALL CORRELATION WITH RST AND DST

- Digit span test measures short-term memory
- Reading span test measures working memory
- Recall scores correlate moderately with both DST and RST, suggesting that Recall measures some of the variances shown in both measures
- Supports use of Recall scores as a screener/ profiler for memory issues



# TOLERABLE TIME WITH NOISE ACCEPTANCE (TNT)

- Tolerable time (willingness to spend time in noise) correlates moderately with TNT (how much noise one is willing to accept without decreasing conversational speech understanding in noise to below 90%).
- TNT has been shown to correlate with hearing aid satisfaction in loud, noisy situations (Seper et al 2019)
- Would TT predict hearing aid satisfaction also?



# CORRELATIONS AMONG VARIOUS OUTCOME MEASURES ON Q-RRT

RRT Outcome	RRT Noise	Recall	Effort	T.Time
Repeat	2TBN-0°	-0.25	-0.66**	-0.42
	SSN-180°	0.22	-0.57*	-0.30
Recall	2TBN-0°	--	0.28	-0.19
	SSN-180°	--	-0.31	0.13
Effort	2TBN-0°	--	--	-0.55*
	SSN-180°	--	--	-0.52*

\*  $p < 0.05$ , \*\*  $p < 0.01$ , \*\*\*  $p < 0.001$ ; Benjamini & Hochberg adjusted.

Repeat scores correlate with listening effort, suggesting LE was based on listeners' listening difficulty

Listening effort correlates with Tolerable time, suggesting more effortful listening, less willingness to spend time

Lack of correlation between Recall and LE/TT/Repeat, suggesting "memory" does not play major role on these tasks (on Q-RRT), and Recall is a separate measure

Repeat does not correlate with TT, suggesting that objective speech understanding alone is not what determines TT

# WHEN SHOULD I USE THE QUICK RRT?

As  
p  
eva  
**UNAI**


After  
**VAI**  
m



ccess


e “blue”  
reen”

# EVALUATING SPEECH IN NOISE NEEDS



- ≥ 95 %-tile
- 75 - < 95 %-tile
- 25 - < 75 %-tile
- 5 - < 25 %-tile
- < 5 %-tile
- Excellent
- Very good
- Good
- OK, improvement desirable
- Poor


# GATHERING INSIGHTS – WORKING MEMORY CAPACITY




- ≥ 95 %-tile
- 75 - < 95 %-tile
- 25 - < 75 %-tile
- 5 - < 25 %-tile
- < 5 %-tile
  
- Excellent
- Very good
- Good
- OK, improvement desirable
- Poor



# GATHERING INSIGHTS – EFFORTFUL LISTENING



# DIFFERENTIATING PATIENTS – WILLINGNESS TO STAY IN NOISE



Patient S21 is easier to manage/more satisfied




- Better WMC/ recall/ cognition
- LE decreases as SNR improves
- Greater willingness to stay in noise (acceptance of noise)

■  $\geq 95$  %-tile  
■ 75 - < 95 %-tile  
■ 25 - < 75 %-tile  
■ 5 - < 25 %-tile  
■ < 5 %-tile


■ Excellent  
■ Very good  
■ Good  
■ OK, improvement desirable  
■ Poor

# MEASURING SUCCESS: UNAIDED VS AIDED

## (QUICK RRT)






# COUNSELING: AS LEAD-IN TO RECOMMENDATIONS (UNAIDED OR AIDED OWN AID)



# IMPROVING SPEECH IN NOISE

- Ensuring a proper fit
- Choose higher level technology
- Choose ones with speech enhancement algorithms
  - Directional microphones more at poorer SNRs
  - Speech enhancer optimizes comfort at SNR around 10 dB while ensuring speech intelligibility index
  - Also wind noise and impulse noise








# PEOPLE WITH BETTER vs POORER RECALL

- **Better** repeat in noise, lower listening effort, longer tolerable time
- Both groups benefit from use of **HD Locator**
- People with **poorer recall** reported **less effort** with the **Speech Enhancer**
- Speech enhancement features are a MUST for people with poor working memory
- Planning for areas for **rehabilitation**

# COUNSELING: MOTIVATION, DEMONSTRATION AND SETTING REALISTIC EXPECTATIONS (AIDED MOMENT)




# LEAD-IN TO INSTRUCTIONS ON FEATURE: DIRECTIONAL FOCUS





# SUMMARY: REPEAT-RECALL TEST

- **Measures more than speech-in-noise**
- **Efficient** (6 min for Quick-RRT; 20 minutes for full RRT)
- **Characterization** (or profiling) of patients
  - Pinpoint areas of strength & weakness
  - Direct focus for rehabilitation
  - Grouping patients for research
- **Validation** - Demonstration of hearing aid (feature) benefits
- **Counseling** - Comparison with norms to set realistic expectations



# THE REPEAT- RECALL TEST (RRT) TUTORIALS

[www.orca-us.info](http://www.orca-us.info)

- Rationale and development
- Operations
- Clinical applications
- Research applications

Thank you for listening!!!

Visit our website ([www.orca-us.info](http://www.orca-us.info))  
Email: [Francis.kuk@wsa.com](mailto:Francis.kuk@wsa.com) for  
link to software/collaborations

